



Similar picture


## HEWI active+ Grab bar 477.31D100

- rail with right-angled ends and fixing roses
- provides support and extra safety in the bathtub
- for wall mounting with roses
- fixing with roses
- 250 mm wide (centre-to-centre), 85 mm deep and 28 mm in diameter, 60 mm rose diameter
- made of high-quality polyamide in HEWI colour 98 (signal white), with active antimicrobial microsilver



HEWI *active* <sup>+</sup>

## Dimensions in mm



Colour

98

Technical information subject to alteration, 16th August 2019

HEWI Heinrich Wilke GmbH  
Postfach 1260  
D-34442 Bad Arolsen

Telefon: +49 5691 82-366  
Telefax: +49 5691 82-99366  
international@hewi.com  
www.hewi.com

HEWI (UK) Limited  
Holm Oak Barn, Beluncle Halt  
Stoke Road  
Hoo Kent ME3 9NT

Phone +44 1634 258 200  
Fax +44 1634 250 099  
info@hewi.co.uk  
www.hewi.com